# N95 Extreme Comfort™ Particulate Respirator

All products on this page are

MADE  
IN U.S.A.

## 3230 - N95 RESPIRATOR

Up to 50% lower breathing resistance compared to many molded N95 designs.

NIOSH N95 provides at least 95% filtration efficiency against certain non-oil based particles.

### 3230 N95 NIOSH Approved TC-84A-9285

#### EXTREME COMFORT™

- Extremely low breathing resistance of 4.1 mm H<sub>2</sub>O inhalation and 5.1 mm exhalation - you hardly notice you're wearing an N95 respirator
- Latex-free synthetic rubber head straps for comfortable tension profile
- Staple-free smoothness
- No eyewear fogging
- Soft inner lining

#### HIGHEST QUALITY

- Fluid Resistant – 120 mm HG
- Meets CDC guidelines for TB exposure control
- Individually packaged for hygienic handling

#### TRUSTED

- GERSON® has been manufacturing NIOSH Approved respirators for more than 35 years.
- ISO 9001: 2015 certified manufacturing of respirator protective devices
- Made in Middleboro, MA

Individually  
Wrapped  
Packages

## 3230 - N95 RESPIRATOR

DESCRIPTION	QUANTITY #	MODEL #	PART #
N95 Respirator	50/box, 4 boxes/case	3230	083230

# GERSON®

## Molded vs. Foldable Comparison Chart

### MOLDED

Typical Molded N95 (non-valved)

### FOLDABLE

3230 Series N95 (non-valved)

		✓ NIOSH Approved Various	✓ NIOSH Approved TC-84A-9285
Protection	#3230 NIOSH Approved N95 Respirator	✓ NIOSH Approved Various	✓ NIOSH Approved TC-84A-9285
	#3230+ Approved for Medical Use	X / ✓	TC-84A-9317
Construction	Durable outer shell	✓	✓
	Low profile, provides optimal field of vision	✓	✓
	Staple free construction	X / ✓	✓
	Nosewire enclosed	X	✓
Comfort	Comfort Level	Moderate	Extreme Comfort
	Up to 50% lower breathing resistance compared to many molded N95 designs.	X	✓
	Extremely reduces eyewear fogging	X	✓
Convenience	Folded, pouch style design for easy carrying and compact storage	X	✓
	Individually packaged to help maintain cleanliness and reduce the risk of contamination.	X	✓
	Easy to carry	X	✓
	Vending-Ready	X	✓

## GERSON®

16 Commerce Blvd.

Middleboro, MA 02346

Tel 508.947.4000 Fax 508.947.5442

sales@gersonco.com • www.gersonco.com

© 2021 Louis M. Gerson Co., Inc. All Rights Reserved

# GERSON®

## Molded vs. Foldable Testimonials

### Testimonials

MOLDED	VS.	FOLDABLE
<i>"Molded N95s are what I used to wear. They are effective BUT..."</i>		<i>"I tried Gerson's Foldable N95! This is the mask I will be wearing in 2021 because..."</i>
<i>"I get face fatigue"</i>		<i>"it is extremely comfortable"</i>
<i>"the straps can be tight"</i>		<i>"the straps are great!"</i>
<i>"the straps are a problem"</i>		<i>"the straps have a nice give and take balance"</i>
<i>"it is difficult to wear all day"</i>		<i>"I can wear it all day - no problem!"</i>
<i>"the dome easily crushes - can't take it anywhere"</i>		<i>"before I need to use it, the folded mask fits in my pocket"</i>

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
1. Unwrap respirator and bend nosepiece into a curve to open mask.
2. Separate the pre-scored elastic headband into two individual bands.
3. Hold the headbands with your index fingers and thumbs.
4. Cup the open mask **securely** under your chin.
5. Pull both headbands over your head.
6. Rest the top headband on the crown of your head and bring the bottom one down under your ears.
7. Press the nosepiece around the contour of your nose and cheekbones. Adjust the edges of the mask to make sure there are no gaps.
8. Perform a Fit Check by gently cupping both hands around the mask, then exhale sharply to see if you sense or detect any leaks.

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**WARNING:** A respirator does not eliminate the risk of contracting any disease or infection. Occupational use of this respirator requires that you be fit tested by a qualitative or quantitative fit test and that your employer comply with OSHA 29 CFR 1910.134, which includes medical evaluation, training and fit testing.

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#### WARNINGS

*Read and understand warnings on insert before using.*

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